



MED SUPPLY LIFE

Handheld Ultrasound
Authorized Distributor

DOUBLE PROBE ULTRASOUND

AFFORDABLE
VISIBLE - FAST - EASY

MED SUPPLY LIFE CO., LTD.


Sales | Import | Installation | Maintenance of Medical Imaging Equipment

Model: D5CE - Handheld Convex + Transvaginal Ultrasound

Tel / WhatsApp: +66 84-199-9893

LINE Official: @777dzaoj

E-mail: msupplylife@gmail.com



YOU CAN USE IT ANYTIME & ANYWHERE



APPLICATION

Application Sites:

- OB/GYN
- Urology

Application Directions:

Follicular Monitoring, Prostate/Pelvics, Uterus, Bladder/Urinary Function.



FEATURES

- Composed of two probes, a scanner can complete surface and vaginal examination as your choice.
- Small and light, suitable for bedside operation and monitoring.
- Wireless, by using the disposable protecting cover, it can easily solve the sterilization issue of the probe.



SPECIFICATION

- Scan Angle: Convex 60, Transvaginal 10mm
- Frequency: Convex probe 3.5/5.0 MHz, Transvaginal probe 6.0/8.0 MHz
- Scan Mode: B, B/M, Color, PW
- Display Depth: Convex 90~305mm, Transvaginal 40~100mm
- Power: Built-in lithium battery, wireless charge 2h, work 2.5h
- Working System: iOS, Android, Windows
- WiFi Type: 802.11n / 2.4G / 5G
- Dimension: 350x70x23mm
- Weight: 300g



WATERPROOF

Reduce medical accidents caused by cross infection, protect doctors and patients, and reduce hospital related expenses.



WIRELESS CHARGE

Easy and connect with a wireless power bank.



EASY TO USE

Simple operation, short learning curve for people without ultrasound experience.



WiFi CONNECT

WiFi connect with your Android phone, iPhone or iPad.

HIGH QUALITY IMAGES

See The Difference Great Image Quality Makes

